```
Goal goal = game.getGoal();
if (!goal.isReached(game)) {
  // increment time
 long elapsedTime = game.getElapsedTime() + increment;
  game.setElapsedTime(elapsedTime);
  // equipment variables
 Patient patient = game.getPatient();
  Equipment equipment = game.getEquipment();
 PhysiologicMonitorComponent physiologicMonitor = (PhysiologicMonitorComponent)equipment
       .getComponent(PhysiologicMonitorComponent.class);
  CDIMonitorComponent cdiMonitor = (CDIMonitorComponent)equipment
       .getComponent(CDIMonitorComponent.class);
  PumpComponent pump = (PumpComponent)equipment
       .getComponent(PumpComponent.class);
  TubeComponent tube = (TubeComponent)equipment
       .getComponent(TubeComponent.class);
  OxigenatorComponent oxigenator = (OxigenatorComponent)equipment
       .getComponent(OxigenatorComponent.class);
  VentilatorComponent ventilator = (VentilatorComponent)equipment
       .getComponent(VentilatorComponent.class);
  HeaterComponent heater = (HeaterComponent)equipment
       .getComponent(HeaterComponent.class);
 BubbleDetectorComponent bubbleDetector = (BubbleDetectorComponent)equipment
       .getComponent(BubbleDetectorComponent.class);
  PressureMonitorComponent pressureMonitor = (PressureMonitorComponent)equipment
       .getComponent(PressureMonitorComponent.class);
  AlarmIndicatorComponent alarmIndicator = (AlarmIndicatorComponent)equipment
       .getComponent(AlarmIndicatorComponent.class);
  // update equipment (physiologic monitor)
 physiologicMonitor.setHeartRate(patient.getHeartRate());
 physiologicMonitor.setO2Saturation(patient.getO2Saturation());
 physiologicMonitor.setDiastolicBloodPressure(patient.getDiastolicBloodPressure());
  physiologicMonitor.setSystolicBloodPressure(patient.getSystolicBloodPressure());
  // update equipment (tubing)
  double cdiPo2 = cdiMonitor.getPO2();
  double tubeSa02 = 0;
  if (cdiPo2 != 0) {
     / see (b) in "Oxyhemoglobin Dissociation Curve.xls"
    tubeSaO2 = 1 / ((23400 / ((cdiPo2 * cdiPo2 * cdiPo2) + (150 * cdiPo2))) + 1);
  tube.setSaO2(tubeSaO2);
  tube.setSvO2(tube.getSaO2() * 0.75);
  tube.setPreMembranePressure((pump.getFlow() * 400) + (oxigenator.getClotting() * 50));
  tube.setPostMembranePressure(tube.getPreMembranePressure() - 40);
  tube.setVenousPressure(-100);
  else {
    tube.setVenousPressure(0);
  if (!tube.isArterialAOpen()) {
    tube.setArterialBubbles(true);
  if ((!tube.isArterialBOpen()) && (!tube.isBridgeOpen())) {
    tube.setArterialBubbles(true);
  if ((tube.getVenousPressure() < -75) && (pump.isOn())) {</pre>
    tube.setVenusBubbles(true);
  else if (pump.isOn()) {
   tube.setVenusBubbles(false);
  // update equipment (pressure monitor)
 pressureMonitor.setPreMembranePressure(tube.getPreMembranePressure());
 pressureMonitor.setPostMembranePressure(tube.getPostMembranePressure());
  pressureMonitor.setVenousPressure(tube.getVenousPressure());
  // update equipment (bubble detector)
 bubbleDetector.setAlarm(tube.isArterialBubbles());
  // update equipment (CDI monitor)
 cdiMonitor.setSaO2(tube.getSaO2());
 cdiMonitor.setSvO2(tube.getSvO2());
 cdiMonitor.setTemperature(heater.getTemperature());
 cdiMonitor.setHct(patient.getHct());
 cdiMonitor.setHgb(patient.getHgb());
 cdiMonitor.setHCO3(patient.getHCO3());
```

```
cdiMonitor.setBE(patient.getBE());
cdiMonitor.setPH(patient.getPH());
cdiMonitor.setPO2((oxigenator.getFiO2() * 600) + (oxigenator.getO2Flow() * 10));
cdiMonitor.setPCO2(patient.getPCO2());
 / update equipment (pump)
if (bubbleDetector.isAlarm() || pressureMonitor.isAlarm()) {
  // turn off the pump
  pump.setOn(false);
if ((pump.getFlow() < (0.02f * patient.getWeight())) || (!pump.isOn())) {</pre>
 pump.setAlarm(true);
// update equipment (alarm)
boolean alarm = false;
if (pressureMonitor.isAlarm()) {
  alarm = true;
if (bubbleDetector.isAlarm()) {
  alarm = true;
if (pump.isAlarm()) {
  alarm = true;
if (oxigenator.isAlarm()) {
  alarm = true;
if (heater.isAlarm()) {
  alarm = true;
alarmIndicator.setAlarm(alarm);
double po2 = oxigenator.getFi02() * oxigenator.getO2Flow() * pump.getFlow() * 100;
double patientSaturation = 0;
if (po2 != 0) {
  // see (b) in "Oxyhemoglobin Dissociation Curve.xls"  
  patientSaturation = 1 / ((23400 / ((po2 * po2 * po2) + (150 * po2))) + 1);
patient.set02Saturation(patientSaturation);
patient.setPCO2(40 + oxigenator.getCo2Flow() * 25);
patient.setPO2((oxigenator.getFiO2() * (760 - 47)) - (patient.getPCO2() / 0.8));
```